



# **Integra<sup>®</sup> MicroFrance<sup>®</sup>**

## **Integra ENT - Throat catalogue**

|                                                  |       |                                                     |       |
|--------------------------------------------------|-------|-----------------------------------------------------|-------|
| <b>Table of contents</b>                         | ii    | Tracheal retractor                                  | 16    |
| <b>Warranty and Customer Service Information</b> | vi    | Grooved director                                    | 16    |
| <b>Indications / Contraindications</b>           | vii   | Curved scissors                                     | 16    |
| <b>Diagnosis</b>                                 | 2     | <b>Micro surgery</b>                                | 17-18 |
| Tongue depressor                                 | 2     | Micro irrigation tube                               | 17    |
| <b>Kilner mouth-gag</b>                          | 3     | Micro vascular clamp                                | 17    |
| <b>Davis-Boyle mouth-gag</b>                     | 4     | Micro dissecting forceps                            | 17    |
| <b>Mouth-gag, complete</b>                       | 5-6   | Jeweler forceps                                     | 17    |
| <b>Morinière mouth-gag</b>                       | 7-9   | Approximator                                        | 17    |
| Presentation                                     | 7     | Needle-holder                                       | 18    |
| Composition of the mouth-gag                     | 8     | Micro scissors                                      | 18    |
| Add-ins                                          | 9     | Silicone cover                                      | 18    |
| <b>Tongue depressor</b>                          | 10    | <b>Suction tube with finger control</b>             | 19    |
| Tongue depressor with suction tube and basket    | 10    | Micro suction irrigation tube                       | 19    |
| Tongue depressor with basket                     | 10    | Jankowski                                           | 19    |
| <b>Tonsil adenoid forceps</b>                    | 11    | Lempert                                             | 19    |
| Tonsil forceps with teeth                        | 11    | Wertheimer                                          | 19    |
| Tonsil haemostatic forceps                       | 11    | <b>Elevator - Raspatory</b>                         | 20    |
| Adenoid forceps                                  | 11    | <b>Scalpel handle</b>                               | 21    |
| <b>Tonsil adenoid instruments</b>                | 12-13 | <b>Curved raspatory</b>                             | 22    |
| Scalpel handle                                   | 12    | <b>Other raspatory</b>                              | 23    |
| Tonsil knife                                     | 12    | Sharp raspatory                                     | 23    |
| Blunt curette                                    | 12    | Curved raspatory                                    | 23    |
| Curved elevator                                  | 12    | Straight raspatory                                  | 23    |
| Double-ended elevator                            | 12    | <b>Retractor</b>                                    | 24-25 |
| Suction tube with finger control                 | 12    | Sharp hook                                          | 24    |
| Elevator and dissector                           | 12    | Double retractor                                    | 24    |
| Monopolar elevator suction tube                  | 12    | Blunt retractor                                     | 24    |
| Monopolar electrode pick                         | 12    | Sharp retractor, articulated                        | 24    |
| Monopolar cable                                  | 12    | Retractor with prongs                               | 25    |
| Curved elevator                                  | 13    | <b>Large self-retaining retractor</b>               | 26    |
| <b>Adenoid curette</b>                           | 14    | <b>Articulated self-retaining retractor</b>         | 27    |
| Twist basket                                     | 14    | <b>Malleable blade retractor - Double retractor</b> | 28    |
| without basket                                   | 14    | <b>Blunt retractor - Double retractor</b>           | 29    |
| <b>Tonsil snare, Amygdalotome</b>                | 15    | <b>Small scissors, pointed, matt finish</b>         | 30    |
| Amygdalotome, tonsil blade                       | 15    | <b>Long scissors, pointed</b>                       | 31    |
| Amygdalotome, handle                             | 15    | <b>Small scissors, blunt, matt finish</b>           | 32    |
| Tonsil snare                                     | 15    | <b>Long curved scissors, blunt</b>                  | 33    |
| Snare wire                                       | 15    | <b>Small Metzenbaum scissors, matt finish</b>       | 34    |
| Stainless steel wire                             | 15    | <b>Large Metzenbaum scissors</b>                    | 35    |
| <b>Neck surgery</b>                              | 16    | <b>Mayo scissors, matt finish</b>                   | 36    |
| Facial nerve dissector                           | 16    |                                                     |       |
| Tracheal dilator                                 | 16    |                                                     |       |
| Curved dissector                                 | 16    |                                                     |       |
| Shear                                            | 16    |                                                     |       |

|                                            |       |                                                    |       |
|--------------------------------------------|-------|----------------------------------------------------|-------|
| <b>Fine dissecting forceps</b>             | 37    | Diverticuloscope                                   | 56    |
| with teeth and platform                    | 37    | Suction tube                                       | 56    |
| with serrated jaws                         | 37    | Fiber                                              | 56    |
| with teeth                                 | 37    |                                                    |       |
| <b>Adson dissecting forceps</b>            | 38    | <b>Laryngoscope holder, Cold Light</b>             | 57    |
| <b>Medium dissecting forceps</b>           | 39-40 | Laryngoscope holder                                | 57    |
| with serrated jaws                         | 39-40 | Plate for Chest holder                             | 57    |
| with teeth                                 | 39-40 | Operating table clamp                              | 57    |
| with teeth and platform                    | 39-40 | Cold light cable                                   | 57    |
|                                            |       | Adapter for cold light cable                       | 57    |
| <b>Dissecting forceps, 180mm</b>           | 41    | <b>Open laryngoscope</b>                           | 58    |
| <b>Dissecting forceps, 200mm</b>           | 42    | <b>Laryngoscope without suction tube</b>           | 59    |
| <b>Medium bipolar dissection forceps</b>   | 43    | <b>Laryngoscope with suction tube</b>              | 60    |
| <b>Bipolar dissecting forceps, 200mm</b>   | 44    | <b>Sataloff laryngoscope</b>                       | 61    |
| <b>Irrigation bipolar forceps</b>          | 45    | <b>Micro suction tube</b>                          | 62    |
| <b>Medium monopolar dissecting forceps</b> | 46    | Suction tube                                       | 62    |
| Adson                                      | 46    | Blunt eyes suction tube                            | 62    |
| with smooth jaws                           | 46    | Insulated suction tube with finger control         | 62    |
| with serrated jaws                         | 46    | Suction tube with finger control                   | 62    |
| with teeth                                 | 46    | Sataloff suction tube                              | 62    |
| Monopolar cable                            | 46    |                                                    |       |
| <b>Large monopolar dissecting forceps</b>  | 47    | <b>Bouchayer Instruments</b>                       | 63    |
| <b>Needle-holder, tungsten jaws</b>        | 48    | Forceps, serrated jaws                             | 63    |
| <b>Large tungsten needle-holder</b>        | 49    | Heart-shaped forceps                               | 63    |
| <b>Fine haemostatic forceps</b>            | 50    | Cup forceps                                        | 63    |
| <b>Medium haemostatic forceps</b>          | 51    | <b>Bouchayer Instruments, Working length 220mm</b> | 64    |
| <b>Strong haemostatic forceps</b>          | 52    | Micro laryngeal scissors                           | 64    |
| With teeth, «T» tip                        | 52    | Lancet knife                                       | 64    |
| with serrated jaws                         | 52    | Lancet knife, handle and 3 blades                  | 64    |
| <b>Other haemostatic forceps</b>           | 53    | Micro laryngeal elevator                           | 64    |
| <b>Esophagoscope</b>                       | 54    | Micro retractor forceps                            | 64    |
| Esophagoscope                              | 54    | Needle-holder                                      | 64    |
| Fiber                                      | 54    |                                                    |       |
| Handle                                     | 54    | <b>Pérouse Instruments, Working length 220mm</b>   | 65    |
| Suction tube                               | 54    | Micro laryngeal elevator                           | 65    |
| Forceps                                    | 54    | Stylet handle                                      | 65    |
| Sterilization case                         | 54    | Tungsten needle-holder insert                      | 65    |
| <b>Hypopharyngoscope</b>                   | 55    | <b>Abitbol Instruments</b>                         | 66    |
| Hypopharyngoscope                          | 55    | Insulated suction tube with finger control         | 66    |
| Fiber                                      | 55    | Monopolar forceps with suction                     | 66    |
| Suction tube                               | 55    | Monopolar cable                                    | 66    |
| Forceps                                    | 55    |                                                    |       |
| <b>Diverticuloscope</b>                    | 56    | <b>Other laryngeal instruments</b>                 | 67-68 |
|                                            |       | Scalpel handle                                     | 67    |
|                                            |       | Knife blade                                        | 67    |
|                                            |       | Stylet handle                                      | 67    |
|                                            |       | Tungsten needle-holder insert                      | 67    |
|                                            |       | Adapter for cleaning insert                        | 67    |
|                                            |       | Needle-Holder, tungsten jaws                       | 67    |
|                                            |       | Sickle knife                                       | 67    |

|                                                |       |                                                          |       |
|------------------------------------------------|-------|----------------------------------------------------------|-------|
| Knife                                          | 67    | <b>Henry set - 3.0mm</b>                                 | 75    |
| Hook                                           | 67    | Scissors                                                 | 75    |
| Retractor                                      | 67    | Monopolar hook                                           | 75    |
| Knife                                          | 68    | Palpation probe                                          | 75    |
| Micro scissors                                 | 68    | Forceps, fenestrated jaws                                | 75    |
| Laryngeal sprinkler                            | 68    | Forceps, alligator serrated jaws                         | 75    |
| Injection needle                               | 68    | Trocar                                                   | 75    |
| Double probe                                   | 68    | Trocar holder ring                                       | 75    |
| <b>Micro coagulation</b>                       | 69    | <b>Terris thyroidectomy Instrument Set -</b>             |       |
| Monopolar suction tube with finger control     | 69    | <b>Robotic thyroid surgery</b>                           | 76    |
| Monopolar electrode pick                       | 69    | Double retractor                                         | 76    |
| Monopolar cable                                | 69    | Curved elevator                                          | 76    |
| Bipolar forceps                                | 69    | Peanut holder                                            | 76    |
| Bipolar cable                                  | 69    | Malleable suction tube elevator                          | 76    |
| <b>Micro injection</b>                         | 70    | Robotic thyroid surgery                                  | 76    |
| Casiano transcervical laryngeal injection      | 70    | <b>Mini invasive instrumentation Ø 3.5mm</b>             | 77    |
| Pistol for injection                           | 70    | Metallic trocar with silicone valve                      | 77    |
| Luer-lock adapter                              | 70    | Metallic trocar with magnetic valve                      | 77    |
| Injection needle                               | 70    | <b>Pediatric Instrumentation Ø 3.5mm</b>                 | 78-86 |
| <b>Sataloff instruments</b>                    | 71-72 | Composite trocar with silicone valve                     | 78    |
| Injection needle                               | 71    | Suction irrigation tube                                  | 79    |
| Micro knife                                    | 71    | Retractable knife                                        | 79    |
| Micro blunt hook                               | 71    | Monopolar hook                                           | 79    |
| Micro curette for cyst                         | 71    | Bipolar forceps                                          | 80    |
| Micro spatula                                  | 71    | 2-piece : Insert Jaws with insulated tube                | 81    |
| Micro hook                                     | 71    | Cleaning adapter                                         | 81    |
| Mirror                                         | 71    | Angled handle                                            | 82    |
| Micro ball dissector                           | 71    | Axial handle                                             | 83    |
| Instrument holder                              | 71    | Ergonomic handle                                         | 84    |
| Cotton carrier                                 | 71    | Scissors, 1-piece                                        | 85    |
| Cup forceps                                    | 72    | Forceps, 1-piece                                         | 85    |
| Heart-shaped grasper                           | 72    | De Lagausie instruments                                  | 86    |
| Alligator forceps                              | 72    | Needle-holder                                            | 86    |
| <b>Sataloff indirect Laryngeal instruments</b> | 73    | <b>Stainless steel cup, Tray</b>                         | 87    |
| Injection needle                               | 73    | <b>Stainless steel sterilization case</b>                | 88    |
| Laryngeal dilator                              | 73    | <b>Plastic sterilization case</b>                        | 89    |
| Micro suction tube                             | 73    | <b>Large plastic sterilization case</b>                  | 90    |
| Micro knife                                    | 73    | <b>Numerical index</b>                                   | 91    |
| Micro scissors                                 | 73    | <b>Table for corresponding manufacturer part numbers</b> |       |
| Micro vascular knife                           | 73    | <b>For information purposes only</b>                     | 99    |
| Micro forceps                                  | 73    |                                                          |       |
| Instrument holder                              | 73    |                                                          |       |
| Stainless steel sterilization case             | 73    |                                                          |       |
| <b>Benninger instruments</b>                   | 74    |                                                          |       |
| Scissors                                       | 74    |                                                          |       |
| Cup forceps                                    | 74    |                                                          |       |
| Heart-shaped grasper                           | 74    |                                                          |       |
| Needle-holder                                  | 74    |                                                          |       |
| Knot pusher                                    | 74    |                                                          |       |



## Integra Limited Warranty

Please refer to the warranty included in the packaging of each product to know the applicable warranty.

## Contacts Customer service:

### International

+33 (0) 4 37 47 59 50  
+33 (0) 4 37 47 59 25 (Fax)  
csemea@integralife.com

### France

+33 (0) 4 37 47 59 10  
+33 (0) 4 37 47 59 29 (Fax)  
custservfrance@integralife.com

### Belgium & Luxembourg

+32 (0) 2 257 4130  
+32 (0) 2 253 2466 (Fax)  
custsvcbenelux@integralife.com

### United Kingdom

+44 (0) 1264 312 725  
+44 (0) 1264 312 821 (Fax)  
custsvcs.uk@integralife.com

### Netherlands

+31(0) 852 083 167  
+31(0) 207 093 627 (Fax)  
custsvcnetheerlands@integralife.com

### Italy

+39 (0) 2 577 8921  
+39 (0) 2 575 113 71 (Fax)  
custsvcitaly@integralife.com

## MicroFrance® ENT Instrumentation


### INDICATIONS

- All the medical devices showed in this catalogue have indications for diagnosis and/or ENT surgery.
- Sizes are not contractual and are given for information purposes only.

### CONTRAINDICATIONS

- Do not use in indications other than the indications mentioned in the DFU or Labeling.

## How to read the catalogue



## Explanation of the table reference and complementary product description

| Title or subtitles               |                             |
|----------------------------------|-----------------------------|
| Global complementary description |                             |
| Product ref.1                    | complementary description 1 |
| Product ref.2                    | complementary description 2 |

## Tongue depressor

Doyen



Double

150mm - Pack 10

7798T<sup>(D)\*</sup>

Doyen



Double

150mm - Pack 10

7798T1<sup>(D)\*</sup>

Canuyt



Fenestrated

185mm

7796

<sup>(D)</sup> Distributed by Integra MicroFrance SAS (France) - For corresponding manufacturer part numbers, see table in annex  
\*Please check with your sales representative availability of this product in your country



Scale 2/3



### Kilner

Ring blade

| 7803S-3     | 7803S-2     | 7803S-1     |
|-------------|-------------|-------------|
| 40x75mm n°3 | 35x70mm n°2 | 30x60mm n°1 |



Palatine blade

**7803B**

without fixing upper jaw



### Kilner

Frame

**7803C**



Camo® dental mouth guard

**7960T40<sup>(D)</sup>\***

Pack 40 - Single use



### Chabolle

Ring blade

| MCR3A       | MCR2A       |
|-------------|-------------|
| 40x90mm n°3 | 40x80mm n°2 |



Silicone dental mouth guard

**7802P**

Pack 2

<sup>(D)</sup> Distributed by Integra MicroFrance SAS (France) - For corresponding manufacturer part numbers, see table in annex  
\*Please check with your sales representative availability of this product in your country



### Davis-Boyle

Frame

7802C



Camo® dental mouth guard

7960T40<sup>(D)\*</sup>

Pack 40 - Single use



Silicone dental mouth guard

7802P

Pack 2



### Davis-Boyle


Blade with suction tube

| 7802S-5      | 7802S-4     | 7802S-3     | 7802S-2     | 7802S-1     |
|--------------|-------------|-------------|-------------|-------------|
| 20x105mm n°5 | 25x90mm n°4 | 25x75mm n°3 | 25x60mm n°2 | 20x25mm n°1 |

<sup>(D)</sup> Distributed by Integra MicroFrance SAS (France) - For corresponding manufacturer part numbers, see table in annex

\*Please check with your sales representative availability of this product in your country

Scale 1/2



### Masson

7809MF

140mm



7821B

### Doyen

7821B

120mm

7821A

140mm

Scale 1/2


### Seemann-Seiffert

|         |         |
|---------|---------|
| 7805-11 | 7805-13 |
| 110mm   | 130mm   |





### Jennings


|         |       |
|---------|-------|
| 7804-9  | 90mm  |
| 7804-11 | 110mm |
| 7804-13 | 130mm |



Presentation of Morinière mouth-gag for mini invasive surgery



Scale 1/2



Frame MCL206C  
with ring blade MCL206A1  
and adapter MCL206P

## Composition of the mouth-gag


Scale 1/3

**MCL206**

Mouth-gag, mini invasive surgery, complete with 1 frame, 5 blades, 1 locking key and 1 angled intermediate piece.

Frame

**MCL206C**



Blade

**MCL206A1** 35x120mm ring blade



**MCL206A2** 35x120mm



**MCL206A3** 35x110mm



**MCL206A4** 40x85mm



**MCL206A5** 35x65mm



Locking key

Angled intermediate piece

**MCL206K**



**MCL206P**



Additional accessories, not included in ref. MCL206

**MCL206A6** larynx blade  
140x20mm



**MCL206A7** Zenker's diverticulum  
blade 180x23mm



## Add-ins instruments

Picture for information, without scale

### Dental mouth guard

**7802P** Silicone dental mouth guard (Pk 2)

**7960T40<sup>(D)\*</sup>** Camo® dental mouth guard  
Single use - (Pk 40)

### Fiberoptic

**MCL38N<sup>(D)\*</sup>** Fiberoptic, 2200mm without adapter

### Adapters

**MCL170<sup>(D)(i)\*</sup>**

**MCL171<sup>(D)(i)\*</sup>**

**MCL173<sup>(D)(i)\*</sup>**


**MCL174<sup>(D)(i)\*</sup>**

**MCL197<sup>(D)(i)\*</sup>**

**Wolf**

**Storz**

**Olympus**



<sup>(i)</sup> Not a medical product


### Light fiber

**MCL89-2** Light fiber

### Chest holder

**MCL76** Laryngoscope arm holder without MCL35


**MCL35** Plate for chest holder



### Suction tube

**MCL89-1** Suction tube for MCL206

**MCL100** Blunt eyes suction tube




### Container

**7874MF-34** Stainless steel sterilization case  
380x205x85mm 340x200x80mm




### Silicone cover

**MCM25-3** Silicone cover



<sup>(D)</sup> Distributed by Integra MicroFrance SAS (France) - For corresponding manufacturer part numbers, see table in annex  
\*Please check with your sales representative availability of this product in your country

Portmann



200mm

7857

Lecouteur



160mm

7800MF

Le Mée



7797-2

170mm

7797-1 small

7797-2 large



Scale 1/2

Birkett



Schmidt




Négus




Bourgeois



Portmann



Chatelier



Scale 1/2




**Jost**  
200mm  
7879



**Carpenter**  
220mm  
7815



**Hurd**  
220mm  
7814M



**Hurd**  
220mm  
7814



180mm  
7878

**Monopolar cables**  
3.5m



**CP391E**

Erbe ACC ICC/VIO  
Storz Wolf



**CP391M**

Aesculap Berchtold  
Lamidey  
Martin Erbe T Serie  
Sutter BM-780 II



**CP391V**

Bowa Bovie  
Conmed Valleylab



**Blondiau**  
230mm  
7812  
15x20mm



**Ruault**  
210mm  
7811



**Lérault**  
200mm  
7884  
5mm



220mm  
7895



200mm  
1003LI<sup>(D)\*</sup>  
for blades 10/15

<sup>(D)</sup> Distributed by Integra MicroFrance SAS (France) - For corresponding manufacturer part numbers, see table in annex  
\*Please check with your sales representative availability of this product in your country

## Hurd

Elevator  
225mm

3734042



Scale 2/3

**Moure**

**Moure**

**Lermoyez**



7828-1-13



7828-16




7829-10


220mm - without basket



7828-1-10 6/10mm



7828-1-13 8/13mm



7828-1-16 10/16mm




7828-1-19 12/19mm


220mm - with basket



7828-10 6/10mm



7828-13 8/13mm



7828-16 10/16mm



7828-19 12/19mm

220mm - without basket



7829-10 5/10mm

7829-13 7/13mm


7829-16 9/16mm

7829-19 13/19mm


Scale 1/2

Working length «L»


## Sluder



## Brünings



## Vacher



**Laborde**



Tracheal dilator

7840

120mm 80°

**Guerrier**



Tracheal retractor

7838

180mm

**Marie**



Director

7839

145mm



Facial nerve dissector

7830M

160mm

**Moure**



Shear serrated blades

7835

200mm

**Guerrier**



«S» curved scissors

7921

180mm

**O'Shaugnessy**



Curved dissector

7837-18

180mm

**Redon**



7920

140mm


Working length «L»

Micro irrigation tube




|           |         |        |
|-----------|---------|--------|
| MCM10-0-5 | Ø 0.5mm | L 30mm |
| MCM10-1-0 | Ø 1mm   |        |

Approximator L 13mm




Large

MCM13



Small

MCM13-1



MCM13

Micro dissecting forceps - 150mm




MCM4BISR



MCM4BISR

1x2 teeth and platform



MCM4R

Smooth jaws with plate

Jeweler forceps



130mm

|          |     |
|----------|-----|
| MCM973D3 | n°3 |
| MCM973D5 | n°5 |

Micro scissors

Vanas



**MCM7BISR**

8mm  
Blunt curved

Vanas



**MCM7R**

8mm - serrated  
straight blades



**MCM970**

8mm  
atraumatic curved



**MCM7BISR**  
150mm

Micro needle-holder

O'Brien



**MCM12R**

5x0.3mm  
smooth curved jaws

Barraquer-Comtet




**MCM1R**

8mm  
smooth straight jaws




**MCM12R**  
150mm



**MCM25-3**

210x130mm





- 1.5mm
- 2mm
- 2.5mm
- 3mm
- 3.5mm
- 4mm

|            |       |
|------------|-------|
| MCO201-2-5 | 2.5mm |
| MCO201-3-0 | 3mm   |

|            |       |
|------------|-------|
| Wertheimer |       |
| MCO770-1-5 | 1.5mm |
| MCO770-2-0 | 2mm   |
| MCO770-2-5 | 2.5mm |
| MCO770-3-0 | 3mm   |
| MCO770-3-5 | 3.5mm |
| MCO770-4-0 | 4mm   |

|            |     |
|------------|-----|
| Jankowski  |     |
| MCEN770J40 | 4mm |
| MCEN770J30 | 3mm |

|            |       |
|------------|-------|
| Lempert    |       |
| MCO771-2-5 | 2.5mm |
| MCO771-3-0 | 3mm   |
| MCO771-4-0 | 4mm   |

|           |     |
|-----------|-----|
| CP330-3-0 | 3mm |
| CP330-4-0 | 4mm |

Sébileau



190mm

CP985

5mm

Fleury



180mm

MCO305

8.5mm

Freer



185mm

CP611

5mm

Howarth



220mm

CP625

Cottle



220mm

CP626

McKenty



240mm

CP627



120mm

160mm

200mm

150mm

120mm

1003IG<sup>(D)\*</sup>

1005I<sup>(D)\*</sup>

1003LI<sup>(D)\*</sup>

CP821

1004IG<sup>(D)\*</sup>

N° 3

N° 7

N° 3L

N° 3

N° 4

For blades 11/15

For blades 23

<sup>(D)</sup> Distributed by Integra MicroFrance SAS (France) - For corresponding manufacturer part numbers, see table in annex  
\*Please check with your sales representative availability of this product in your country

Farabeuf



150mm

7787-10

10mm

Fleury



170mm

MCO979

6,5mm

Sébileau



150mm

CP825

4mm

Plester



180mm

MCO988

7,5mm



195mm

3724059

6,5mm

## Joseph




CP310



180mm


CP310

5mm



180mm

CP310D

5mm  
straight

180mm

CP310E


4mm

## Obwegeser



CP344A

## Obwegeser



180mm

CP344A

5mm


CP344B

7mm

CP344C

11mm

## Obwegeser-Tessier



180mm - sharp

CP344AT

5mm

CP344BT

7mm

CP344CT

11mm

## Obwegeser



190mm

CP344D

7mm  
curved

CP344E

7mm  
very curved

## Obwegeser-Tessier



190mm - tranchante

CP344DT

7mm  
curved

CP344ET

7mm  
very curved

## Farabeuf



|        |       |
|--------|-------|
| 140mm  | 120mm |
| 7768-1 | 7768  |
| Pk 2   | Pk 2  |

## Desmarres



|          |
|----------|
| 140mm    |
| CP822-12 |
| 12mm     |
| CP822-10 |
| 10mm     |
| CP822-14 |
| 14mm     |

## Gillies



|        |
|--------|
| 180mm  |
| CP17   |
| 3mm    |
| CP17-1 |
| 2mm    |

## Gillies



|        |
|--------|
| 145mm  |
| CP17-3 |
| 2mm    |
| CP17-2 |
| 3mm    |

## Humby



|             |
|-------------|
| 110mm       |
| CP941       |
| articulated |

## Knapp



150mm

3714027

8mm

3714028

blunt



140mm

CP21-3

pointu

CP21A-3

blunt

## Volkmann



150mm

7782-1-2



180mm

CP335

## Senn-Miller



160mm

CP609

3 sharp prongs

## Weitlaner



curved

160mm

0-28-7

sharp prongs




**Bizeau**

165mm

MCL359

blades 20x30mm



**Weitlaner-Beckmann**

165mm

7781




Malleable



210mm

CP818-8



8mm

CP818-8

10mm

CP818-10

12mm

CP818-12

15mm

CP818-15

20mm

CP818-20

Dautrey



15x40mm

15x55mm

140mm

CP984

Dautrey



15x60mm

15x85mm

140mm

CP984-1

Malleable



10mm

15mm

Blade retractor  
10/15mm - 160mm

7773

Scale 2/3



200mm

CEV418

Pack 2




Descottes

200mm

CEV417

Pack 2

Mosnier



160mm

CEV419

Pack 2

Pessey



180mm

CP335P

Leroux-Robert



160mm

7923

Scale 3/4

## Streiff



CP16-1

105mm

| CP16-1        | CP16            |
|---------------|-----------------|
| curved blades | straight blades |



100mm

| CP353            |
|------------------|
| Strong curvature |

## Offret



100mm

| CP356 |
|-------|
|-------|

## Joseph



CP354A

120mm

| CP354A        | CP354           |
|---------------|-----------------|
| curved blades | straight blades |

## Joseph



CP919A

110mm

| CP919A        | CP919           |
|---------------|-----------------|
| curved blades | straight blades |

## Fox



CP358A

110mm

| CP358A        | CP358           |
|---------------|-----------------|
| curved blades | straight blades |

Kilner-jost



120mm

CP362

Curved

Boetcher



180mm

CP997

Curved

Brophy



CP367-1

150mm

CP367-1

Curved

CP367

Straight

Joseph



CP355A

150mm

CP355A

Curved

CP355

Straight

Scale 3/4

**Shea**



120mm

**CP637**

curved blades

**Bregeat**



110mm

**CP363**

curved blades

**Stevens**



110mm

**CP922A**

curved tenotomy scissors

**Stevens**



110mm

**CP922**

straight blades

**Stevens**




115mm

**CP922B**

curved tenotomy scissors

Reynolds



160mm  
CP921

Ragnell



150mm  
CP616-15



130mm  
CP616-13

Joseph




150mm  
CP355C

Léger




150mm  
CP365



110mm


|               |                 |
|---------------|-----------------|
| CP351-7       | CP351-6         |
| curved blades | straight blades |



150mm

|               |                 |
|---------------|-----------------|
| CP351-3       | CP351-2         |
| curved blades | straight blades |

## Delicate blades



150mm

|               |                 |
|---------------|-----------------|
| CP931-1       | CP931           |
| curved blades | straight blades |





CP351-5

Curved standard

| CP351-5 | CP351-1 |
|---------|---------|
| 200mm   | 180mm   |



CP931-5


Curved, delicate blades

| CP931-5 | CP931-3 |
|---------|---------|
| 200mm   | 180mm   |



## Dissecting forceps 110mm

Castroviejo



CP938-1

Horizontal 1x2 teeth  
and platform

Castroviejo



CP37

Oblique 1x2 teeth  
and platform



CP332

Serrated jaws  
1.4mm

## Dissecting forceps 120mm



serrated jaws

CP414

Jaws 1mm



7831A-12

Jaws 1.5mm



CP406

1x2 teeth - Jaws 1mm



Jaws with teeth

7833A-12

1x2 teeth - Jaws 1.5mm




CP400-1

5x5 griffes

150mm

Adson



|          |                         |
|----------|-------------------------|
| CP604-15 | serrated<br>jaws 1.5mm  |
| CP605-15 | 1x2 teeth<br>1.5mm jaws |

Adson-Brown



|          |           |
|----------|-----------|
| CP606-15 | 7x7 teeth |
|----------|-----------|

120mm



Adson-Tubiana

Jaws 1mm

M832

serrated



M833

1x2 teeth



Adson

Jaws 1.5mm

CP604

serrated



CP605

1x2 teeth



CP606A

tungsten serrated

Adson-Brown

Jaws 7x7 teeth

CP606

## Mac-Indoe

140mm




**CP934**  
Oblique  
1x2 teeth  
and  
serrated  
platform

**CP934-1**  
Serrated  
jaws

## Debakey

## Barraya



150mm

140mm

| CP415                | 7831A-15               | CP407                 | 7833A-15                | CP988-15               | CP933-14                         | CP400     |
|----------------------|------------------------|-----------------------|-------------------------|------------------------|----------------------------------|-----------|
| serrated<br>jaws 1mm | serrated<br>jaws 1.5mm | 1x2 teeth<br>jaws 1mm | 1x2 teeth<br>jaws 1.5mm | serrated<br>jaws 1.5mm | 3x4 teeth &<br>serrated platform | 5x5 teeth |

Précelle



150mm

CP345

serrated  
angled jaws

Barraya



160mm

CP933-16

3x4 teeth and  
serrated platform

Claoué



160mm

CP401

5x6 teeth

Mac-Indoe



170mm

CP936

1x2 teeth and  
serrated platform




## Debakey

## Barraya



1.5mm  
7831A-18



1.5mm  
CP988-18



1mm  
CP409  
1x2



1.5mm  
7833A-18  
1x2



1.5mm  
CP933-18  
3x4 teeth and  
serrated platform



6x6  
CP400-2

Debaquey



| 200mm           | 210mm            |
|-----------------|------------------|
| 7832            | 7834             |
| 1mm<br>serrated | 1mm<br>1x2 teeth |



|          |
|----------|
| 2mm      |
| 7831A-20 |



|          |
|----------|
| 1,5mm    |
| CP988-20 |



|           |
|-----------|
| 2mm       |
| 7833A-20  |
| 1x2 teeth |



Bipolar cables - 4.5 m



CP393-2  
Bowa Emc Lamidey Valleylab



CP393-3  
Storz Wolf Erbe (except ICC/VIO international)



CP393-4  
Aesculap Berchtold Martin  
Sutter BM-780 II

Magnan



Adson



| 110mm                |                        | 120mm              |                      | 160mm                |                        | 200mm                  |  |
|----------------------|------------------------|--------------------|----------------------|----------------------|------------------------|------------------------|--|
| CP392-1R             | CP604BR                | CP605BR            | CP392-2R             | CP392-8R             | CP392BI                | CP392BS                |  |
| smooth jaws<br>1x8mm | serrated Jaws<br>1.5mm | 1x2 teeth<br>1.5mm | smooth jaws<br>1x8mm | smooth jaws<br>1x8mm | smooth jaws<br>0.5x5mm | smooth jaws<br>0.5x5mm |  |



Bipolar cables - 4.5 m



**CP393-2**

Bowa Emc Lamidey Valleylab



**CP393-3**

Storz Wolf Erbe  
(except ICC/VIO  
international)



**CP393-4**

Aesculap Berchtold Martin  
Sutter BM-780 II

**CP392-3R**

Straight  
smooth jaws  
1x8mm

**CP392-10R**

Angled tip  
smooth jaws  
1x8mm

**CP392-5R**

Bayonet  
smooth jaws  
1x8mm

**CP392-9R**

Bayonet  
Angled down jaws  
1x8mm

Dousset



300mm

CEV392D30

Jaws 3mm

Belghiti



230mm

CEV392D20

Jaws 2mm



230mm

CEV392P20

Jaws 2mm



Soubrane



230mm

CEV392S20

Jaws 2mm

Gayet



240mm

CEV392G24

Jaws 5x15mm

Partenski



210mm

CEV392R15

Jaws 4x10mm

Bipolar cables - 4.5 m



CP393-2  
Bowa Emc Lamidey Valleylab



CP393-3  
Storz Wolf Erbe  
(except ICC/VIO international)



CP393-4  
Aesculap Berchtold Martin  
Sutter BM-780 II

Scale 2/3

Working length - «L»

## Chambon



290mm L 260mm

**CP390L26**

smooth jaws 1x8mm



**CP390-3**

smooth jaws 1x8mm



230mm L 200mm

**CP390-5**

angled jaws 1x8mm



**CP390-7**

smooth jaws 1x8mm

Monopolar cables 3.5m



**CP391E**

Erbe ACC ICC/VIO  
Storz Wolf



**CP391M**

Aesculap Berchtold Martin  
Erbe T Serie Lamidey  
Sutter BM-780 II



**CP391V**

Bowa Bovie Conmed Valleylab

Monopolar cables - 3.5 m



**CP391E**  
Erbe ACC ICC/VIO  
Storz Wolf



**CP391M**  
Aesculap Berchtold Lamidey  
Martin Erbe T Serie  
Sutter BM-780 II



**CP391V**  
Bowa Bovie  
Conmed Valleylab

## Adson



110mm

**CP390-1**

smooth jaws 1x8mm



120mm

**CP604M**

serrated jaws 1.5mm



**CP605M**

1x2 teeth 1.5mm



160mm

**CP390-2**

smooth jaws 1x8mm




**CP390-8**

1x2 teeth 1x8mm



### Halsey



130mm

| 703T          | 704T        |
|---------------|-------------|
| serrated jaws | smooth jaws |



703T

### Crile-Wood Mayo-Hégar




150mm

| 702T | 714T |
|------|------|
|------|------|




702T

### Becker




130mm

| CP1107       |
|--------------|
| diamond jaws |



Scale 2/3


Crile-Wood



200mm  
720T



Crile-Wood




180mm  
701T



Stratte



210mm  
7945T



## Halsted

120mm - with teeth

|        |        |
|--------|--------|
| P761-7 | P761-8 |
|--------|--------|

|          |        |
|----------|--------|
| Straight | Curved |
|----------|--------|



P761-7

## Leriche

140mm - without teeth

|        |        |
|--------|--------|
| 7867-1 | 7867-2 |
|--------|--------|

|          |        |
|----------|--------|
| Straight | Curved |
|----------|--------|

140mm - with teeth

|        |        |
|--------|--------|
| 7867-3 | 7867-4 |
|--------|--------|

|          |        |
|----------|--------|
| Straight | Curved |
|----------|--------|



7867-1



Kelly

without teeth

| 140mm  | 160mm  |
|--------|--------|
| P762-2 | P762-4 |
| Curved | Curved |



P762-2

Kocher


140mm 160mm 180mm

Straight without teeth

|        |        |   |
|--------|--------|---|
| P764-1 | P764-5 | - |
|--------|--------|---|

Straight with teeth

|        |        |         |
|--------|--------|---------|
| P764-3 | P764-7 | P764-11 |
|--------|--------|---------|



P764-3



### Martel

130mm

P763

Forceps with teeth  
'T' tip



### Terrier

130mm

P769

Haemostatic forceps  
serrated jaws



### Péan

140mm 160mm

P765-14

P765-16

Haemostatic forceps  
serrated jaws

P765-14




### Chaput

Forceps with teeth

130mm

7866-13

2x3 teeth



Haemostatic  
fenestrated forceps  
160mm

P770

Heart serrated jaws




### Ombredanne

with teeth  
130mm

P771

## Blunt eyes suction tube

Scale 1/4



Working length - «L»

MCL940L45

|         | 2mm       | 2.5mm     | 3mm       | 3.5mm     | 4mm       |
|---------|-----------|-----------|-----------|-----------|-----------|
| L 450mm | MCL920L45 | MCL925L45 | MCL930L45 | MCL935L45 | MCL940L45 |
| L 550mm | MCL920L55 | MCL925L55 | MCL930L55 | MCL935L55 | MCL940L55 |

## Forceps





7844T5L40



| Forceps | Crocodile forceps | Cup forceps<br>Ø 4mm | Cup forceps up opening<br>3.5x10mm | Cup forceps upward<br>Ø 4mm |
|---------|-------------------|----------------------|------------------------------------|-----------------------------|
| L 400mm | 7844T3L40         | 7844T5L40            | 7844T6L40                          | 7844T7L40                   |
| L 450mm | 7844T3L45         | 7844T5L45            | 7844T6L45                          | 7844T7L45                   |
| L 500mm | 7844T3L50         | 7844T5L50            | 7844T6L50                          | 7844T7L50                   |
| L 550mm | 7844T3L55         | 7844T5L55            | 7844T6L55                          | 7844T7L55                   |
| L 600mm | 7844T3L60         | 7844T5L60            | 7844T6L60                          | 7844T7L60                   |

## Esophagoscope

**MCL320** Esophagoscope  
without handle - fiber



**MCL330** Esophagoscope  
without handle - fiber



Sterilization case

**7874MF-61**

640x189x95mm - L 610x180x80mm

Scale 1/3

Scale 1/4

Working length - «L»





## Blunt eyes suction tube



|         | 2mm       | 2.5mm     | 3mm       | 3.5mm     | 4mm       |           |
|---------|-----------|-----------|-----------|-----------|-----------|-----------|
| L 250mm | MCL920L25 | MCL925L25 | MCL930L25 | MCL935L25 | MCL940L25 | MCL930L25 |
| L 350mm | MCL920L35 | MCL925L35 | MCL930L35 | MCL935L35 | MCL940L35 |           |


## Forceps

| Forceps | Crocodile forceps | Cup forceps<br>Ø 4mm | Cup forceps up opening<br>3.5x10mm | Cup forceps upward<br>Ø 4mm |           |
|---------|-------------------|----------------------|------------------------------------|-----------------------------|-----------|
| L 200mm | 7844T3L20         | 7844T5L20            | 7844T6L20                          | 7844T74L20                  | 7844T5L25 |
| L 250mm | 7844T3L25         | 7844T5L25            | 7844T6L25                          | 7844T74L25                  |           |
| L 300mm | 7844T3L30         | 7844T5L30            | 7844T6L30                          | 7844T74L30                  |           |
| L 350mm | 7844T3L35         | 7844T5L35            | 7844T6L35                          | 7844T74L35                  |           |




## Hypopharyngoscope







## Lacau diverticuloscope



## Reyt diverticuloscope



Picture for information, without scale



<sup>(1)</sup> Not a medical product

Cold light





| Adapters       | generator side                                                                                                | instrument side                                                                                                 |
|----------------|---------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------|
| <b>Wolf</b>    | MCL170 <sup>(D)(1)*</sup>  | MCL171 <sup>(D)(1)*</sup>  |
| <b>Storz</b>   | MCL173 <sup>(D)(1)*</sup>  | MCL174 <sup>(D)(1)*</sup>  |
| <b>Olympus</b> | MCL197 <sup>(D)(1)*</sup>  |                                                                                                                 |

<sup>(D)</sup> Distributed by Integra MicroFrance SAS (France) - For corresponding manufacturer part numbers, see table in annex  
\*Please check with your sales representative availability of this product in your country


Scale 1/3

Working length - «L»


## Lescanne




## Garabédian




## Garcin



## Legent



## Legent






Scale 1/3

Working length - «L»

## Bouchayer

**MCL79**  
without fiber




L 160mm

Fiber

**MCL78**

## Bouchayer


**MCL77**  
without fiber



L 160mm

## Pignat

**MCL64**  
without fiber  
small size




L 120mm

Fiber

**MCL64F**

## Pignat

**MCL63**  
without fiber  
large size





L 180mm


Fiber

**MCL63F**


Scale 1/3





Abitbol



Brasnu



Abitbol




Scale 1/3

Sataloff

Working length - «L»



**MCL120**

Pediatric size  
With fiber and suction





**MCL121**

Medium male size  
With fiber and suction





**MCL122**

Adult size  
With fiber and suction





**MCL123**

Medium female size  
With fiber and suction





**MCL124**

«Vallecula»  
With fiber and suction



**MCL125**

With fiber and suction



| Complete laryngoscope | Right light fiber | Left light fiber | Dual suction tube |
|-----------------------|-------------------|------------------|-------------------|
| MCL120                | MCLS120-2         | MCLS120-1        | MCLS120-3         |
| MCL121                | MCLS121-2         | MCLS121-1        | MCLS121-3         |
| MCL122                | MCLS122-2         | MCLS122-1        | MCLS122-3         |
| MCL123                | MCLS123-2         | MCLS123-1        | MCLS1233          |
| MCL124                | MCL124FR          | MCL124FL         | MCL124ST          |
| MCL125                | MCL125FR          | MCL125FL         | MCL125ST          |

## Suction tube

**Bouchayer**



MCL10

Ø 2.5mm

L 220mm

Working length - «L»

Scale 1/3

## Blunt eyes suction tube

**Brasnu**



MCL100

Ø 2.5mm

L 220mm

**Garabédian**



MCL100T20

Ø 2mm


MCL100T15

Ø 1.5mm

L 220mm

## Insulated suction tube

**Abitbol**



MCL210

Ø 2.5mm

L 220mm

Stop and go

## Dual suction tube

**Huttenbrink**



MCL103

Ø 2.5mm

L 240mm

Stop and go

## French velvet eye suction tube

**Sataloff**



Micro suction tube

MCL-S16



Micro suction tube  
with retractor below

MCL-S17



Micro suction tube with  
retractor above

MCL-S17-1

300mm - Ø 1mm


L 230mm

Stop and go

MCL-S17

Working length - «L»

Scale 1/3



**MCL18**

curved right

**MCL15**


curved left

**MCL25**

straight



## Alligator forceps




**MCL19**


left

**MCL19BIS**

right



## Heart-shaped forceps 45°



**MCL12**

straight

**MCL29**


curved left  
45°

**MCL30**

curved right  
45°

**MCL33**

upward



## Cup forceps Ø 2mm



**MCL67**


upward

**MCL68**

curved left

**MCL69**

curved right




## Oval cup forceps


Scale 1/2

## Micro scissors

### Bouchayer



## Lancet knife, elevator



## Micro needle-holder




## Micro retractor forceps



Scale 1/2


**Elevator**



## Stylet handle and needle-holder insert

Scale 1/3

Micro needle-holder  
handle and insert



**MCEN39L**


Luer-lock adapter for  
tube insert cleaning port  
Pack 3

Stylet handle without ratchet

**MCL828**




## Insulated suction tube with finger control



MCL210

Ø 2.5mm



## Monopolar fenestrated forceps



MCL205

3x9mm jaws



## Heart-shaped forceps, monopolar



MCL219

right curved



MCL219-1

left curved



## Alligator forceps, monopolar




MCL241

right curved



MCL241-1

left curved



Monopolar cables - 3.5 m



CP391E

Erbe ACC ICC/VIO Storz Wolf



CP391M

Aesculap Berchtold Lamidey Martin Erbe T Serie Sutter BM-780 II



CP391V

Bowa Bovie Conmed Valleylab



Scale 1/2 Working length - «L»

Knife blade sterile, Single use, Pack 25  
Scale 2/3



SM63<sup>(D)\*</sup>



SM64<sup>(D)\*</sup>



SM65<sup>(D)\*</sup>



Laryngeal diastema - pediatrics

CP821E

Scalpel handle 170mm

Scalpel handle, knife blade

Garabédian

Scale 1/3

L 180mm

Micro pediatric needle-holder  
handle and insert

Tungsten needle-holder insert



MCL94IG

delivered with adapter  
MCEN39L for cleaning



MCEN39L

Luer-lock adapter for  
tube insert cleaning port  
Pack 3

Stylet handle with ratchet

MCL828C

Scale 1/2

L 230mm

Needle-holder, tungsten jaws



MCL94G

1-piece

Cleaning  
instrument

Garabédian

Oval knife, hook, retractor

Scale 1/2

L 220mm



MCL4

oval knife 3x4mm



MCL5

oval knife 45° 2mm



MCL7

blunt hook



MCL48

blunt retractor 7mm

Micro sickle knife

Scale 1/2

Micro sickle knife curved 10mm

L 220mm




MCL1BIS



<sup>(D)</sup> Distributed by Integra MicroFrance SAS (France) - For corresponding manufacturer part numbers, see table in annex  
\*Please check with your sales representative availability of this product in your country

## Parsons

Laryngeal sprinkler



**MCL98**

Ø 1,5mm

L 200mm

Scale 1/2  
Working length - «L»

## Bouchayer

Bayonet injection needle



**MCL58**

Ø 0,5mm

L 210mm

## Woo

Scissors, 5mm blades



**MCLW10**


straight



**MCLW11**

up 45°

L 200mm



**3734115**

knife

L 250mm



## Lowinger


Double probe

**MCL101**

Ø 2,5mm / 1mm  
280mm

**MCL102**

Ø 2mm / 1mm  
260mm




Scale 1/3

Working length - «L»

## Frèche

Electrode pick bayonet




MCL59

L 220mm



## Monopolar instruments


Suction tube with finger control



MCL62

Ø 2.5mm

L 220mm



Monopolar cables - 3.5 m



CP391E

Erbe ACC ICC/VIO Storz Wolf



CP391M

Aesculap Berchtold Lamidey Martin Erbe T Serie Sutter BM-780 II



CP391V

Bowa Bovie Conmed Valleylab

Micro forceps, straight jaws




MCL34



Delivered with box and tray

L 280mm



## Bipolar forceps

## Peretti


Micro forceps, angled jaws



MCL340

Jaws 45° - 2mm

L 210mm



Bipolar cables - 4.5 m



CP393-2

Bowa Emc Lamidey Valleylab



CP393-3

Storz Wolf Erbe (except ICC/VIO international)



CP393-4

Aesculap Berchtold Martin Sutter BM-780 II

## Injector, injection needles

### Bouchayer

Injection needle  
bayonet tip Ø 1mm



**MCL55-1**



**MCL55-2**

right



**MCL55-3**

left

Injection needle, straight  
malleable Ø 0.5mm




**MCL55U**

Pk 6 - Single use  
non sterile



L 220mm



**MCL55-1**

L 245mm



**MCL55U**

Scale 1/3

Working length - «L»

Injector

**MCL55**

with needle MCL55-1

**MCL55B**

with adapter MCL55R

Luer-lock adapter



**MCL55R**

Scale 1/1


## Transcervical laryngeal injection

### Casiano

**MCL55T6**

**MCL55T4**

**MCL55T5**



Handle for injection needle

**MCL55T4**

Trocar obturator


**MCL55T6**

Injection needle  
Ø 0.8mm L 40mm - Pack 6


**MCL55T5**

non sterile

Scale 1/3




## Dismantled instruments for instrument holder MCL-S19



| Instruments L 220mm                                                                 |                                                                                     |                                                                                     |                                                                                     |                                                                                     | L 250mm                                                                               |  |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|--|
| Hook                                                                                | Mirror                                                                              | Ball dissector                                                                      | Cotton carrier                                                                      | Curved curette for cyst                                                             | Mirror                                                                                |  |
|  |  |  |  |  |  |  |
| <b>MCL-S2</b>                                                                       | <b>MCL-S14</b>                                                                      | <b>MCL-S10</b>                                                                      | <b>MCL-S3</b>                                                                       | <b>MCL-S38</b>                                                                      | <b>MCL-S15</b>                                                                        |  |
| 2 prongs, sharp                                                                     | 70° - Ø 3mm                                                                         |                                                                                     | 19mm                                                                                | Ø 1mm                                                                               | 70° - Ø 4mm                                                                           |  |

Scale 2/3  
Working length - «L»



MCL-S18

Cup forceps



**MCL-S18**  
downbiting


Heart-shaped grasper



**MCL-S21**  
right



**MCL-S22**  
left



MCL-S24

Alligator forceps



**MCL-S23**  
right



**MCL-S24**  
left




**MCL-S39**  
fine right jaws



**MCL-S41**  
fine left jaws

Scale 1/3  
Working length - «L»



Cup forceps Ø 2.4mm - L 180mm



MCLS49

straight



MCLS50

curved left



MCLS51

curved right

L 180mm

Heart-shaped forceps 45° L 180mm



MCLS55

Curved left



MCLS56

Curved right

MCLS55

Scissors 3mm blades - L 180mm



MCLS52

straight



MCLS53

curved right



MCLS54

curved left

Alligator forceps Ø 5mm - L 180mm



MCLS57

straight



MCLS58

curved right



MCLS59

curved left

MCLS65

Micro knives 70mm angled - Pack 3  
for instrument holder MCLS48



MCLS63

vascular  
2mm



MCLS64

straight sickle  
7mm



MCLS65

curved sickle  
6mm



MCLS66

Flap  
90° 3mm

Instrument holder

MCLS48



Stainless steel sterilization case

MCLS62SC

L 190mm

Suction tube  
L 180mm



MCLS60

Ø 1.8mm



MCLS61

Ø 2.5mm

MCLS60

Injection needle Ø 0.7mm - Pack 6



MCL-S42

curved



MCL-S43

straight

Scale 1/3

Working length 220mm

Micro scissors



**MCLB16**

curved right




**MCLB17**

curved left



**MCLB31**

angled 70°



**MCLB31**


Cup forceps



**MCLS18R**

280mm

downbiting



Heart-shaped grasper



**MCL-S21R**

right



**MCL-S22R**

left



**MCL-S23R**

right



**MCL-S24R**

left



**MCL-S22R**

Alligator forceps

Needle-driver



**MCLB94**

280mm

angled up 70°



Suture knot pusher



**MCLB30**

250mm







Scale 1/2  
Working length «L»



Trocara holder ring  
Ø 10mm  
**MCLP70**



Trocara  
with blunt tip  
**MCLP10**



Fenestrated forceps, 3x8mm jaws  
**MCLP20U**




Fenestrated forceps, 3x8mm jaws  
**MCLP25U**



forceps, alligator jaws, 6mm  
**MCLP30U**



Micro scissors, 4mm blades  
**MCLP40U**



Monopolar hook  
**MCLP50**



**MCLP50**




Round palpation probe  
**MCLP60**



**MCLP60**

Scale 1/2

|                                                                                     |                                                                                      |
|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| Blunt retractor 180mm                                                               |                                                                                      |
| MCLT10<br>small tips : 8x35 / 15x40mm                                               |    |
| MCLT11<br>Large tips : 8x35 / 15x50mm                                               |                                                                                      |
| Blunt curved elevator 190mm                                                         |                                                                                      |
| MCLT20<br>2mm                                                                       |    |
| MCLT21<br>4mm                                                                       |                                                                                      |
| Suction tube elevator 200mm Ø 3mm                                                   |                                                                                      |
| MCLT40<br>Malleable                                                                 |    |
| Nerve eye suction tube with finger control 200mm Ø 3mm                              |                                                                                      |
| MCLT50<br>Malleable                                                                 |   |
| Peanut holder 200mm Ø 3mm axial handle                                              |                                                                                      |
| MCLT30                                                                              |  |
|  |                                                                                      |


## Robotic thyroid surgery

|                |                                          |                                                                                                       |
|----------------|------------------------------------------|-------------------------------------------------------------------------------------------------------|
| <b>MCLT12R</b> | Muscle retractor                         | <br><br>not shown |
| <b>MCLT13R</b> | Bivalve retractor                        |                                                                                                       |
| <b>MCLT21R</b> | Ball hook elevator Ø 4mm 440mm           |                                                                                                       |
| <b>MCLT40R</b> | Suction elevator, malleable, Ø 3mm 270mm |                                                                                                       |
| <b>MCLT50R</b> | Malleable ner eye suction tube, 330mm    |                                                                                                       |

not shown

## Trocar with metal cannula

|                           |                  |                                                                                      | Working length                                                                  |          |
|---------------------------|------------------|--------------------------------------------------------------------------------------|---------------------------------------------------------------------------------|----------|
|                           |                  |                                                                                      | 50 mm                                                                           | 80 mm    |
| Trocar tip                | smooth           |     | CEV8143A4                                                                       | CEV143A4 |
|                           | pyramidal        |    | CEV8143A3                                                                       | CEV143A3 |
| Sealing cap               | green            |     | CEV103P6<br>Pack 6                                                              |          |
| Trocar head               | magnetic         |    | CEV177A3<br>Delivered complete with 1 magnetic valve<br>and 1 green sealing cap |          |
|                           | silicone         |     | CEV127T3<br>Delivered complete with 1 silicone valve<br>and 1 green sealing cap |          |
| Valve for trocar head     | magnetic         |   | CEV179A5<br>Pack 6                                                              |          |
|                           | silicone         |   | CEV125P6<br>Pack 6                                                              |          |
| Metal cannula, crenelated | with stopcock    |  | CEV8143A7                                                                       | CEV143A7 |
|                           | without stopcock |   | CEV8143B7                                                                       | CEV143B7 |



Scale 1/1


## Composite trocar with silicone valve




|                                                                                              |                                                                                                                          |  | Working length     |             |
|----------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------|--|--------------------|-------------|
|                                                                                              |                                                                                                                          |  | 45 mm              | 70 mm       |
| Trocar tip  | smooth                                  |  | CEV8826M (1)       | CEV326M (1) |
|                                                                                              | pyramidal                               |  | CEV8826P (1)       | CEV326P (1) |
| Sealing cap                                                                                  | green                                   |  | CEV103P6<br>Pack 6 |             |
| Valve for trocar head                                                                        | silicone                                |  | CEV125P6<br>Pack 6 |             |
| Trocar body composed of<br>1 trocar head<br>with composite cannula<br>and silicone valve     | smooth cannula<br>without stopcock     |  | CEV8826L (1)       | CEV326L (1) |
|                                                                                              | threaded cannula<br>without stopcock  |  | CEV8826V (1)       | CEV326V (1) |

(1) New reference


## Monopolar hook




## Retractable knife



## Suction irrigation tube



(1) New reference

| 1-piece                                                                                                       |                                                                                                            |     |                    |
|---------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------|-----|--------------------|
| Complete instruments : bipolar forceps                                                                        |                                                                                                            |     |                    |
| Jaws<br>Scale 1/1                                                                                             | Handle                                                                                                     |     | Working length     |
| <br>serrated jaws, 1.5mm     | <br>Angled, bipolar       |     | CEV897<br>250 mm   |
|                             |                                                                                                            | 1/3 |                    |
| <br>fenestrated jaws, 10x3mm | <br>«K» handle, bipolar   |     | CEV897BG<br>250 mm |
|                           |                                                                                                            | 1/3 |                    |
| <br>spatula jaws, 2mm      | <br>«K» handle, bipolar |     | CEV897C<br>250 mm  |
|                           |                                                                                                            | 1/3 |                    |

| Bipolar cables – 4.5 m                                                              |                                                                                     |                                                                                       |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
|  |  |  |
| CP393-2                                                                             | CP393-3                                                                             | CP393-4                                                                               |
| Bowa Emc Lamidey Valleylab                                                          | Storz Wolf Erbe (except ICC/VIO international)                                      | Aesculap Berchtold Martin Sutter BM-780 II                                            |

| Semi-Modular                                                                                                                                                                                                                               |                                            |                                                                                                                                                                                                                                                                                                                    |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Inserts with tubes                                                                                                                                                                                                                         |                                            | Handle                                                                                                                                                                                                                                                                                                             |
| Jaws with insulated tube<br>Scale 1/1                                                                                                                                                                                                      | Working length                             | Pages                                                                                                                                                                                                                                                                                                              |
| <br>Babcock insert                                                                                                                                        | CEV891IU<br>235 mm                         | <br><br><br><br><br><br><b>Rotatable</b><br>page 82-83 |
| <br>Curved scissors insert serrated blades 8mm                                                                                                            | CEV893IU<br>230 mm                         |                                                                                                                                                                                                                                                                                                                    |
| <br>Straight scissors insert serrated blades 8mm                                                                                                          | CEV894IU<br>230 mm                         |                                                                                                                                                                                                                                                                                                                    |
| <br>Fenestrated insert 3x12mm jaws                                                                                                                       | CEV895IU<br>235 mm                         |                                                                                                                                                                                                                                                                                                                    |
| <br>Grasping insert                                                                                                                                     | CEV899IU<br>220 mm<br><i>New reference</i> |                                                                                                                                                                                                                                                                                                                    |
| <br>Alligator jaws insert 6mm                                                                                                                           | CEV896IU<br>235 mm                         |                                                                                                                                                                                                                                                                                                                    |
| <br>Dissector insert<br>medium curved 10mm jaws                                                                                                         | CEV898IU<br>230 mm                         |                                                                                                                                                                                                                                                                                                                    |
|  <div>CEV39L<br/>Adapter cleaning luer-lock for tube insert</div>  |                                            |                                                                                                                                                                                                                                                                                                                    |

## Modular

### Rotatable instruments

Angled handles with black rotatable ring  
Scale 1/3



with ratchet

CEV1019T5BU  
Large pattern



with ratchet

CEV10195CU



with ratchet

CEV10195DU



with ratchet

CEV10195RU



without ratchet

CEV1039T5TBU  
Large pattern



without ratchet

CEV10395DU



without ratchet

CEV10395RU



without ratchet


CEV1039G5U



Modular

Rotatable instruments


Axial handles with black rotatable ring  
Scale 1/3



CEV6285R



CEV6385RU



CEV628P5



CEV728T5

## Modular

### Rotatable instruments

Ergonomic handles with black rotatable ring  
Scale 1/3



CEV720BRUC

Large ergonomic handle  
with ratchet



CEV720RU

Fine ergonomic handle  
with ratchet



CEV7285RU

**Elhadad**


Fine ergonomic handle  
with ratchet



CEV728BR5U

Large ergonomic handle  
with ratchet

| 1-piece                                                                                                                                       |                                                                                                                 |                                           |
|-----------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------|-------------------------------------------|
| Complete instruments                                                                                                                          |                                                                                                                 |                                           |
| Jaws<br>Scale 1/1                                                                                                                             | Handle                                                                                                          | Working length                            |
| <br><b>Montupet</b><br>Babcock forceps                       | <br>Angled with long ratchet   | <b>CEV891MU</b><br>200 mm                 |
| <br><b>Montupet</b><br>Straight scissors serrated blades 8mm | <br>Angled                     | <b>CEV894MU</b><br>220 mm                 |
| <br><b>Montupet</b><br>Fenestrated forceps 3x12mm jaws      | <br>Angled with long ratchet  | <b>CEV895MU</b><br>195 mm                 |
| <br>Grasping forceps with teeth                            | <br>Angled with long ratchet | <b>CEV899U</b><br>220 mm<br>New reference |

| 1-piece                                                                                                                                          |                                                                                    |       |                           |
|--------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------|-------|---------------------------|
| De Lagausic instruments                                                                                                                          |                                                                                    |       |                           |
| Jaws<br>Scale 1/1                                                                                                                                | Handle                                                                             |       | Working length            |
| <br><b>De Lagausic</b><br>Fenestrated forceps 3x12mm jaws       |   | Axial | <b>CEV995MU</b><br>220 mm |
| <br><b>De Lagausic</b><br>Straight scissors serrated blades 8mm |   | Axial | <b>CEV996MU</b><br>220 mm |
| <br><b>De Lagausic</b><br>Alligator forceps 6mm serrated jaws   |   | Axial | <b>CEV997MU</b><br>210 mm |
| <br><b>De Lagausic</b><br>Dissector medium curved 10mm          |  | Axial | <b>CEV998MU</b><br>220 mm |

| Needle-holder                                                                                                                                  |                                                                                     |                             |                          |
|------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-----------------------------|--------------------------|
| <br><b>Montupet</b><br>Needle-holder straight tungsten jaws |  | Angled<br>with long ratchet | <b>CEV892</b><br>220 mm  |
| <br><b>Reinberg</b><br>Needle-holder straight tungsten jaws |  | Axial                       | <b>CEV992</b><br>220 mm  |
| <br>Needle-holder left curved tungsten jaws                 |  | Axial                       | <b>CEV992D</b><br>220 mm |

Cup



7876-60<sup>(D)\*</sup>

Cup

|                          |       |
|--------------------------|-------|
| 7876-60 <sup>(D)*</sup>  | 60mm  |
| 7876-80 <sup>(D)*</sup>  | 80mm  |
| 7876-100 <sup>(D)*</sup> | 100mm |
| 7876-120 <sup>(D)*</sup> | 120mm |
| 7876-160 <sup>(D)*</sup> | 160mm |

Stainless steel tray




7872-24<sup>(D)\*</sup>  
Scale 1/2

Stainless steel tray

|                         |           |
|-------------------------|-----------|
| 7872-24 <sup>(D)*</sup> | 240x180mm |
|-------------------------|-----------|


<sup>(D)</sup> Distributed by Integra MicroFrance SAS (France) - For corresponding manufacturer part numbers, see table in annex  
\*Please check with your sales representative availability of this product in your country

Scale 1/5




**7874MF-22**

220x100x40mm




**7874MF-25**

250x150x60mm




**7874MF-34**

340x200x80mm



**7874MF-40**

400x200x90mm



**7874MF-61**

610x180x80mm

Picture for information, without scale



MCEN136SC

270x160x35mm



3717011

528x250x123mm




Picture for information, without scale



3717001

530x250x140mm



3717023

525x250x180mm



**Numerical index**

**O**

|              |    |
|--------------|----|
| O-28-7 ..... | 26 |
|--------------|----|

**1**

|                 |        |
|-----------------|--------|
| 1003IG(D) ..... | 21     |
| 1003LI(D) ..... | 12, 21 |
| 1004IG(D) ..... | 21     |
| 1005I(D) .....  | 21     |

**3**

|               |    |
|---------------|----|
| 3714027 ..... | 25 |
| 3714028 ..... | 25 |
| 3717001 ..... | 89 |
| 3717011 ..... | 13 |
| 3717023 ..... | 90 |
| 3724059 ..... | 22 |
| 3734042 ..... | 13 |
| 3734115 ..... | 68 |

**7**

|                  |         |
|------------------|---------|
| 701T .....       | 49      |
| 702T .....       | 48      |
| 703T .....       | 48      |
| 704T .....       | 48      |
| 714T .....       | 48      |
| 720T .....       | 49      |
| 7755-1-0-4 ..... | 15      |
| 7755-1-0-5 ..... | 15      |
| 7755-1-0-6 ..... | 15      |
| 7768 .....       | 24      |
| 7768-1 .....     | 24      |
| 7773 .....       | 28      |
| 7781 .....       | 27      |
| 7782-1-2 .....   | 25      |
| 7787-10 .....    | 22      |
| 7796 .....       | 2       |
| 7797-1 .....     | 10      |
| 7797-2 .....     | 10      |
| 7798T1(D) .....  | 2       |
| 7798T(D) .....   | 2       |
| 7800MF .....     | 10      |
| 7802C .....      | 4       |
| 7802P .....      | 3, 4, 9 |
| 7802S-1 .....    | 4       |
| 7802S-2 .....    | 4       |
| 7802S-3 .....    | 4       |
| 7802S-4 .....    | 4       |
| 7802S-5 .....    | 4       |
| 7803B .....      | 4       |

|                 |    |
|-----------------|----|
| 7803C .....     | 3  |
| 7803S-1 .....   | 3  |
| 7803S-2 .....   | 3  |
| 7803S-3 .....   | 3  |
| 7804-9 .....    | 6  |
| 7804-11 .....   | 6  |
| 7804-13 .....   | 6  |
| 7805-11 .....   | 6  |
| 7805-13 .....   | 6  |
| 7808L-1 .....   | 15 |
| 7808L-2 .....   | 15 |
| 7808L-3 .....   | 15 |
| 7808M .....     | 15 |
| 7809MF .....    | 5  |
| 7810 .....      | 15 |
| 7811 .....      | 12 |
| 7812 .....      | 12 |
| 7813 .....      | 15 |
| 7813-1 .....    | 15 |
| 7814 .....      | 12 |
| 7814M .....     | 12 |
| 7815 .....      | 12 |
| 7819 .....      | 11 |
| 7820 .....      | 11 |
| 7821A .....     | 5  |
| 7821B .....     | 5  |
| 7822 .....      | 11 |
| 7823 .....      | 11 |
| 7825 .....      | 11 |
| 7826-1-7 .....  | 11 |
| 7828-1-10 ..... | 14 |
| 7828-1-13 ..... | 14 |
| 7828-1-16 ..... | 14 |
| 7828-1-19 ..... | 14 |
| 7828-10 .....   | 14 |
| 7828-13 .....   | 14 |
| 7828-16 .....   | 14 |
| 7828-19 .....   | 14 |
| 7829-10 .....   | 14 |
| 7829-13 .....   | 14 |
| 7829-16 .....   | 14 |
| 7829-19 .....   | 14 |
| 7830M .....     | 16 |
| 7831A-12 .....  | 37 |
| 7831A-15 .....  | 39 |
| 7831A-18 .....  | 41 |
| 7831A-20 .....  | 42 |
| 7832 .....      | 42 |
| 7833A-12 .....  | 37 |
| 7833A-15 .....  | 39 |
| 7833A-18 .....  | 41 |
| 7833A-20 .....  | 42 |
| 7834 .....      | 42 |

|                  |       |
|------------------|-------|
| 7835 .....       | 16    |
| 7837-18 .....    | 16    |
| 7838 .....       | 16    |
| 7839 .....       | 16    |
| 7840 .....       | 16    |
| 7843-1 .....     | 55    |
| 7843A-1 .....    | 55    |
| 7844T3L20 .....  | 55    |
| 7844T3L25 .....  | 55    |
| 7844T3L30 .....  | 55    |
| 7844T3L35 .....  | 55    |
| 7844T3L40 .....  | 54    |
| 7844T3L45 .....  | 54    |
| 7844T3L50 .....  | 54    |
| 7844T3L55 .....  | 54    |
| 7844T3L60 .....  | 54    |
| 7844T5L20 .....  | 55    |
| 7844T5L25 .....  | 55    |
| 7844T5L30 .....  | 55    |
| 7844T5L35 .....  | 55    |
| 7844T5L40 .....  | 54    |
| 7844T5L45 .....  | 54    |
| 7844T5L50 .....  | 54    |
| 7844T5L55 .....  | 54    |
| 7844T5L60 .....  | 54    |
| 7844T6L20 .....  | 55    |
| 7844T6L25 .....  | 55    |
| 7844T6L30 .....  | 55    |
| 7844T6L35 .....  | 55    |
| 7844T6L40 .....  | 54    |
| 7844T6L45 .....  | 54    |
| 7844T6L50 .....  | 54    |
| 7844T6L55 .....  | 54    |
| 7844T6L60 .....  | 54    |
| 7844T74L20 ..... | 55    |
| 7844T74L25 ..... | 55    |
| 7844T74L30 ..... | 55    |
| 7844T74L35 ..... | 55    |
| 7844T74L40 ..... | 54    |
| 7844T74L45 ..... | 54    |
| 7844T74L50 ..... | 54    |
| 7844T74L55 ..... | 54    |
| 7844T74L60 ..... | 54    |
| 7857 .....       | 10    |
| 7866-13 .....    | 53    |
| 7867-1 .....     | 50    |
| 7867-2 .....     | 50    |
| 7867-3 .....     | 50    |
| 7867-4 .....     | 50    |
| 7872-24(D) ..... | 87    |
| 7874MF-22 .....  | 88    |
| 7874MF-25 .....  | 88    |
| 7874MF-34 .....  | 9, 88 |

<sup>(D)</sup> Distributed by Integra MicroFrance SAS (France) - For corresponding manufacturer part numbers, see table in annex

|                  |         |
|------------------|---------|
| 7874MF-40.....   | 88      |
| 7874MF-61.....   | 54, 88  |
| 7876-60(D).....  | 87      |
| 7876-80(D).....  | 87      |
| 7876-100(D)..... | 87      |
| 7876-120(D)..... | 87      |
| 7876-160(D)..... | 87      |
| 7878.....        | 12      |
| 7879.....        | 12      |
| 7884.....        | 12      |
| 7895.....        | 12      |
| 7920.....        | 16      |
| 7921.....        | 16      |
| 7923.....        | 29      |
| 7945T.....       | 49      |
| 7960T40(D).....  | 3, 5, 9 |

## CEV

|                 |        |
|-----------------|--------|
| CEV39L.....     | 81     |
| CEV103BG.....   | 79     |
| CEV103GH.....   | 79     |
| CEV103P6.....   | 77, 78 |
| CEV125P6.....   | 77, 78 |
| CEV127T3.....   | 77     |
| CEV143A3.....   | 77     |
| CEV143A4.....   | 77     |
| CEV143A7.....   | 77     |
| CEV143B7.....   | 77     |
| CEV177A3.....   | 77     |
| CEV179A5.....   | 77     |
| CEV326L.....    | 78     |
| CEV326M.....    | 78     |
| CEV326P.....    | 78     |
| CEV326V.....    | 78     |
| CEV392D20.....  | 45     |
| CEV392D30.....  | 45     |
| CEV392G24.....  | 45     |
| CEV392P20.....  | 45     |
| CEV392R15.....  | 45     |
| CEV392S20.....  | 45     |
| CEV417.....     | 29     |
| CEV418.....     | 29     |
| CEV419.....     | 29     |
| CEV628P5.....   | 83     |
| CEV720BRUC..... | 84     |
| CEV720RU.....   | 84     |
| CEV728BR5U..... | 84     |
| CEV728T5.....   | 83     |
| CEV890.....     | 79     |
| CEV890C.....    | 79     |
| CEV891MU.....   | 85     |
| CEV892.....     | 86     |

|                   |    |
|-------------------|----|
| CEV894MU.....     | 85 |
| CEV895MU.....     | 85 |
| CEV896.....       | 79 |
| CEV897.....       | 80 |
| CEV897BG.....     | 80 |
| CEV897C.....      | 80 |
| CEV898.....       | 79 |
| CEV898L.....      | 79 |
| CEV899U.....      | 85 |
| CEV984.....       | 57 |
| CEV992.....       | 86 |
| CEV992D.....      | 86 |
| CEV995MU.....     | 86 |
| CEV996MU.....     | 86 |
| CEV997MU.....     | 86 |
| CEV998MU.....     | 86 |
| CEV1019T5BU.....  | 82 |
| CEV1039G5U.....   | 82 |
| CEV1039T5TBU..... | 82 |
| CEV6285R.....     | 83 |
| CEV6385RU.....    | 83 |
| CEV7285RU.....    | 84 |
| CEV8103G3.....    | 79 |
| CEV8143A3.....    | 77 |
| CEV8143A4.....    | 77 |
| CEV8143A7.....    | 77 |
| CEV8143B7.....    | 77 |
| CEV8826L.....     | 78 |
| CEV8826M.....     | 78 |
| CEV8826P.....     | 78 |
| CEV8826V.....     | 78 |
| CEV10195CU.....   | 82 |
| CEV10195DU.....   | 82 |
| CEV10195RU.....   | 82 |
| CEV10395DU.....   | 82 |
| CEV10395RU.....   | 82 |

## CP

|                |    |
|----------------|----|
| CP16.....      | 30 |
| CP16-1.....    | 30 |
| CP17.....      | 24 |
| CP17-1.....    | 24 |
| CP17-2.....    | 24 |
| CP17-3.....    | 24 |
| CP21-3.....    | 25 |
| CP21A-3.....   | 25 |
| CP37.....      | 37 |
| CP310.....     | 23 |
| CP310D.....    | 23 |
| CP310E.....    | 23 |
| CP330-3-0..... | 19 |
| CP330-4-0..... | 19 |

|               |                    |
|---------------|--------------------|
| CP332.....    | 37                 |
| CP335.....    | 25                 |
| CP335P.....   | 29                 |
| CP344A.....   | 23                 |
| CP344AT.....  | 23                 |
| CP344B.....   | 23                 |
| CP344BT.....  | 23                 |
| CP344C.....   | 23                 |
| CP344CT.....  | 23                 |
| CP344D.....   | 23                 |
| CP344DT.....  | 23                 |
| CP344E.....   | 23                 |
| CP344ET.....  | 23                 |
| CP345.....    | 40                 |
| CP351-1.....  | 35                 |
| CP351-2.....  | 34                 |
| CP351-3.....  | 34                 |
| CP351-5.....  | 35                 |
| CP351-6.....  | 34                 |
| CP351-7.....  | 34                 |
| CP352.....    | 36                 |
| CP352-1.....  | 36                 |
| CP352-2.....  | 36                 |
| CP352-3.....  | 36                 |
| CP353.....    | 30                 |
| CP354.....    | 30                 |
| CP354A.....   | 30                 |
| CP355.....    | 31                 |
| CP355A.....   | 31                 |
| CP355C.....   | 33                 |
| CP356.....    | 30                 |
| CP358.....    | 30                 |
| CP358A.....   | 30                 |
| CP362.....    | 31                 |
| CP363.....    | 32                 |
| CP365.....    | 33                 |
| CP367.....    | 31                 |
| CP367-1.....  | 31                 |
| CP390-1.....  | 47                 |
| CP390-2.....  | 47                 |
| CP390-3.....  | 46                 |
| CP390-5.....  | 46                 |
| CP390-7.....  | 46                 |
| CP390-8.....  | 47                 |
| CP390L26..... | 46                 |
| CP391E.....   | 12, 46, 47, 66, 69 |
| CP391M.....   | 12, 46, 47, 66, 69 |
| CP391V.....   | 12, 46, 47, 66, 69 |
| CP392-1R..... | 43                 |
| CP392-2R..... | 43                 |
| CP392-3R..... | 44                 |
| CP392-5R..... | 44                 |
| CP392-8R..... | 43                 |

<sup>(D)</sup> Distributed by Integra MicroFrance SAS (France) - For corresponding manufacturer part numbers, see table in annex

|                 |                    |
|-----------------|--------------------|
| CP392-9R.....   | 44                 |
| CP392-10R ..... | 44                 |
| CP392BI.....    | 43                 |
| CP392BS.....    | 43                 |
| CP393-2.....    | 43, 44, 45, 69, 80 |
| CP393-3.....    | 43, 44, 45, 69, 80 |
| CP393-4.....    | 43, 44, 45, 69, 80 |
| CP400.....      | 39                 |
| CP400-1 .....   | 37                 |
| CP400-2.....    | 41                 |
| CP401 .....     | 40                 |
| CP406.....      | 37                 |
| CP407 .....     | 39                 |
| CP409.....      | 41                 |
| CP414.....      | 37                 |
| CP415.....      | 39                 |
| CP604.....      | 38                 |
| CP604-15 .....  | 38                 |
| CP604BR.....    | 43                 |
| CP604M.....     | 47                 |
| CP605.....      | 38                 |
| CP605-15.....   | 38                 |
| CP605BR.....    | 43                 |
| CP605M.....     | 47                 |
| CP606 .....     | 38                 |
| CP606-15 .....  | 38                 |
| CP606A.....     | 38                 |
| CP609 .....     | 25                 |
| CP611 .....     | 20                 |
| CP616-13.....   | 33                 |
| CP616-15.....   | 33                 |
| CP625 .....     | 20                 |
| CP626 .....     | 20                 |
| CP627 .....     | 20                 |
| CP628 .....     | 24                 |
| CP637.....      | 32                 |
| CP818-8.....    | 28                 |
| CP818-10 .....  | 28                 |
| CP818-12.....   | 28                 |
| CP818-15.....   | 28                 |
| CP818-20.....   | 28                 |
| CP821.....      | 21                 |
| CP821E.....     | 67                 |
| CP822-10.....   | 24                 |
| CP822-12 .....  | 24                 |
| CP822-14 .....  | 24                 |
| CP825.....      | 22                 |
| CP919 .....     | 30                 |
| CP919A.....     | 30                 |
| CP921.....      | 33                 |
| CP922 .....     | 32                 |
| CP922A.....     | 32                 |
| CP922B.....     | 32                 |
| CP931.....      | 34                 |

|                |    |
|----------------|----|
| CP931-1.....   | 34 |
| CP931-3 .....  | 35 |
| CP931-5 .....  | 35 |
| CP933-14 ..... | 39 |
| CP933-16 ..... | 40 |
| CP933-18 ..... | 41 |
| CP934 .....    | 39 |
| CP934-1 .....  | 39 |
| CP936 .....    | 40 |
| CP938-1 .....  | 37 |
| CP941.....     | 24 |
| CP984 .....    | 28 |
| CP984-1.....   | 28 |
| CP985 .....    | 20 |
| CP988-15.....  | 39 |
| CP988-18.....  | 41 |
| CP988-20 ..... | 42 |
| CP997 .....    | 31 |
| CP1107 .....   | 48 |

## M

|           |    |
|-----------|----|
| M832..... | 38 |
| M833..... | 38 |

## MCEN

|                  |        |
|------------------|--------|
| MCEN39L.....     | 65, 67 |
| MCEN136SC .....  | 89     |
| MCEN770J30 ..... | 19     |
| MCEN770J40.....  | 19     |

## MCL

|               |    |
|---------------|----|
| MCL1BIS.....  | 67 |
| MCL4.....     | 67 |
| MCL5 .....    | 67 |
| MCL7 .....    | 67 |
| MCL10 .....   | 62 |
| MCL11 .....   | 64 |
| MCL12 .....   | 63 |
| MCL15.....    | 63 |
| MCL16 .....   | 64 |
| MCL17.....    | 64 |
| MCL18 .....   | 63 |
| MCL19 .....   | 63 |
| MCL19BIS..... | 63 |
| MCL25 .....   | 63 |
| MCL29 .....   | 63 |
| MCL30.....    | 63 |
| MCL31.....    | 64 |
| MCL33 .....   | 63 |
| MCL34.....    | 69 |

|                |       |
|----------------|-------|
| MCL35 .....    | 9, 57 |
| MCL37 .....    | 58    |
| MCL37C.....    | 58    |
| MCL37CF.....   | 58    |
| MCL37D .....   | 58    |
| MCL37DF .....  | 58    |
| MCL37G .....   | 58    |
| MCL38N(D)..... | 9, 57 |
| MCL42 .....    | 58    |
| MCL46.....     | 60    |
| MCL46A .....   | 60    |
| MCL47 .....    | 60    |
| MCL48.....     | 67    |
| MCL49.....     | 56    |
| MCL49F.....    | 56    |
| MCL51 .....    | 58    |
| MCL52 .....    | 58    |
| MCL53 .....    | 58    |
| MCL54 .....    | 58    |
| MCL55 .....    | 70    |
| MCL55-1.....   | 70    |
| MCL55-2 .....  | 70    |
| MCL55-3 .....  | 70    |
| MCL55B.....    | 70    |
| MCL55R.....    | 70    |
| MCL55T4 .....  | 70    |
| MCL55T5.....   | 70    |
| MCL55T6 .....  | 70    |
| MCL55U .....   | 70    |
| MCL56BL .....  | 64    |
| MCL56BR .....  | 64    |
| MCL56PL .....  | 65    |
| MCL56PR .....  | 65    |
| MCL57 .....    | 64    |
| MCL58 .....    | 68    |
| MCL59 .....    | 69    |
| MCL62 .....    | 69    |
| MCL63.....     | 59    |
| MCL63F.....    | 59    |
| MCL64.....     | 59    |
| MCL64F.....    | 59    |
| MCL65A .....   | 60    |
| MCL65AF .....  | 60    |
| MCL65AS .....  | 60    |
| MCL67 .....    | 63    |
| MCL68.....     | 63    |
| MCL69.....     | 63    |
| MCL76 .....    | 9, 57 |
| MCL77 .....    | 59    |
| MCL78 .....    | 59    |
| MCL79.....     | 59    |
| MCL81BL .....  | 64    |
| MCL81BR.....   | 64    |
| MCL81PL .....  | 65    |

|                   |           |
|-------------------|-----------|
| MCL81PR.....      | 65        |
| MCL82.....        | 64        |
| MCL83.....        | 60        |
| MCL83-1.....      | 60        |
| MCL84.....        | 60        |
| MCL88.....        | 56        |
| MCL88-1.....      | 56        |
| MCL88-2.....      | 56        |
| MCL89.....        | 60        |
| MCL89-1.....      | 9, 56, 60 |
| MCL89-2.....      | 9, 60     |
| MCL93.....        | 64        |
| MCL93-2.....      | 64        |
| MCL94.....        | 64        |
| MCL94G.....       | 67        |
| MCL94IG.....      | 67        |
| MCL94IP.....      | 65        |
| MCL94IPR.....     | 65        |
| MCL98.....        | 68        |
| MCL100.....       | 9, 62     |
| MCL100T15.....    | 62        |
| MCL100T20.....    | 62        |
| MCL101.....       | 68        |
| MCL102.....       | 68        |
| MCL103.....       | 62        |
| MCL110.....       | 55        |
| MCL111.....       | 55        |
| MCL112.....       | 55        |
| MCL112F.....      | 55        |
| MCL112H.....      | 54, 55    |
| MCL120.....       | 61        |
| MCL121.....       | 61        |
| MCL122.....       | 61        |
| MCL123.....       | 61        |
| MCL124.....       | 61        |
| MCL124FL.....     | 61        |
| MCL124FR.....     | 61        |
| MCL124ST.....     | 61        |
| MCL125.....       | 61        |
| MCL125FL.....     | 61        |
| MCL125FR.....     | 61        |
| MCL125ST.....     | 61        |
| MCL170(D)(1)..... | 9, 57     |
| MCL171(D)(1)..... | 9, 57     |
| MCL173(D)(1)..... | 9, 57     |
| MCL174(D)(1)..... | 9, 57     |
| MCL197(D)(1)..... | 9, 57     |
| MCL205.....       | 66        |
| MCL206.....       | 8         |
| MCL206A1.....     | 7, 8      |
| MCL206A2.....     | 8         |
| MCL206A3.....     | 8         |
| MCL206A4.....     | 8         |
| MCL206A5.....     | 8         |

|                |        |
|----------------|--------|
| MCL206A6.....  | 8      |
| MCL206A7.....  | 8      |
| MCL206C.....   | 7, 8   |
| MCL206K.....   | 8      |
| MCL206P.....   | 7, 8   |
| MCL210.....    | 62, 66 |
| MCL219.....    | 66     |
| MCL219-1.....  | 66     |
| MCL241.....    | 66     |
| MCL241-1.....  | 66     |
| MCL290.....    | 60     |
| MCL291.....    | 60     |
| MCL291AD.....  | 60     |
| MCL291AG.....  | 60     |
| MCL291FD.....  | 60     |
| MCL291FG.....  | 60     |
| MCL294.....    | 60     |
| MCL295.....    | 60     |
| MCL320.....    | 54     |
| MCL320F.....   | 54     |
| MCL330.....    | 54     |
| MCL330F.....   | 54     |
| MCL340.....    | 69     |
| MCL359.....    | 27     |
| MCL828.....    | 65     |
| MCL828C.....   | 67     |
| MCL920L25..... | 55     |
| MCL920L30..... | 56     |
| MCL920L35..... | 55     |
| MCL920L45..... | 54     |
| MCL920L55..... | 54     |
| MCL925L25..... | 55     |
| MCL925L30..... | 56     |
| MCL925L35..... | 55     |
| MCL925L45..... | 54     |
| MCL925L55..... | 54     |
| MCL930L25..... | 55     |
| MCL930L30..... | 56     |
| MCL930L35..... | 55     |
| MCL930L45..... | 54     |
| MCL930L55..... | 54     |
| MCL935L25..... | 55     |
| MCL935L30..... | 56     |
| MCL935L35..... | 55     |
| MCL935L45..... | 54     |
| MCL935L55..... | 54     |
| MCL940L25..... | 55     |
| MCL940L30..... | 56     |
| MCL940L35..... | 55     |
| MCL940L45..... | 54     |
| MCL940L55..... | 54     |

## MCLB

|             |    |
|-------------|----|
| MCLB16..... | 74 |
| MCLB17..... | 74 |
| MCLB30..... | 74 |
| MCLB31..... | 74 |
| MCLB94..... | 74 |

## MCLP

|              |    |
|--------------|----|
| MCLP10.....  | 75 |
| MCLP20U..... | 75 |
| MCLP25U..... | 75 |
| MCLP30U..... | 75 |
| MCLP40U..... | 75 |
| MCLP50.....  | 75 |
| MCLP60.....  | 75 |
| MCLP70.....  | 75 |

## MCL-S

|                |    |
|----------------|----|
| MCL-S0.....    | 71 |
| MCL-S1.....    | 71 |
| MCL-S2.....    | 71 |
| MCL-S3.....    | 71 |
| MCL-S4.....    | 71 |
| MCL-S5.....    | 71 |
| MCL-S6.....    | 71 |
| MCL-S7.....    | 71 |
| MCL-S8.....    | 71 |
| MCL-S9.....    | 71 |
| MCL-S10.....   | 71 |
| MCL-S11.....   | 71 |
| MCL-S12.....   | 71 |
| MCL-S13.....   | 71 |
| MCL-S14.....   | 71 |
| MCL-S15.....   | 71 |
| MCL-S16.....   | 62 |
| MCL-S17.....   | 62 |
| MCL-S17-1..... | 62 |
| MCL-S18.....   | 72 |
| MCL-S18R.....  | 74 |
| MCL-S19.....   | 71 |
| MCL-S20.....   | 71 |
| MCL-S21.....   | 72 |
| MCL-S21R.....  | 74 |
| MCL-S22.....   | 72 |
| MCL-S22R.....  | 74 |
| MCL-S23.....   | 72 |
| MCL-S23R.....  | 74 |
| MCL-S24.....   | 72 |
| MCL-S24R.....  | 74 |
| MCL-S25.....   | 71 |
| MCL-S26.....   | 71 |

<sup>(1)</sup> Not a medical product

<sup>(1D)</sup> Distributed by Integra MicroFrance SAS (France) - For corresponding manufacturer part numbers, see table in annex

|                |    |
|----------------|----|
| MCL-S27.....   | 71 |
| MCL-S35.....   | 71 |
| MCL-S38.....   | 71 |
| MCL-S39.....   | 72 |
| MCL-S41.....   | 72 |
| MCL-S42.....   | 73 |
| MCL-S43.....   | 73 |
| MCLS44.....    | 73 |
| MCLS45.....    | 73 |
| MCLS46.....    | 73 |
| MCLS47.....    | 73 |
| MCLS48.....    | 73 |
| MCLS49.....    | 73 |
| MCLS50.....    | 73 |
| MCLS51.....    | 73 |
| MCLS52.....    | 73 |
| MCLS53.....    | 73 |
| MCLS54.....    | 73 |
| MCLS55.....    | 73 |
| MCLS56.....    | 73 |
| MCLS57.....    | 73 |
| MCLS58.....    | 73 |
| MCLS59.....    | 73 |
| MCLS60.....    | 73 |
| MCLS61.....    | 73 |
| MCLS62SC.....  | 73 |
| MCLS63.....    | 73 |
| MCLS64.....    | 73 |
| MCLS65.....    | 73 |
| MCLS66.....    | 73 |
| MCLS120-1..... | 61 |
| MCLS120-2..... | 61 |
| MCLS120-3..... | 61 |
| MCLS121-1..... | 61 |
| MCLS121-2..... | 61 |
| MCLS121-3..... | 61 |
| MCLS122-1..... | 61 |
| MCLS122-2..... | 61 |
| MCLS122-3..... | 61 |
| MCLS123-1..... | 61 |
| MCLS123-2..... | 61 |
| MCLS1233.....  | 61 |

## MCLT

|              |    |
|--------------|----|
| MCLT10.....  | 76 |
| MCLT11.....  | 76 |
| MCLT12R..... | 76 |
| MCLT13R..... | 76 |
| MCLT20.....  | 76 |
| MCLT21.....  | 76 |
| MCLT21R..... | 76 |

|              |    |
|--------------|----|
| MCLT30.....  | 76 |
| MCLT40.....  | 76 |
| MCLT40R..... | 76 |
| MCLT50.....  | 76 |
| MCLT50R..... | 76 |

## MCLW

|             |    |
|-------------|----|
| MCLW10..... | 68 |
| MCLW11..... | 68 |

## MCM

|                |       |
|----------------|-------|
| MCM1R.....     | 18    |
| MCM4BISR.....  | 17    |
| MCM4R.....     | 17    |
| MCM7BISR.....  | 18    |
| MCM7R.....     | 18    |
| MCM10-0-5..... | 17    |
| MCM10-1-0..... | 17    |
| MCM12R.....    | 18    |
| MCM13.....     | 17    |
| MCM13-1.....   | 17    |
| MCM25-3.....   | 9, 18 |
| MCM970.....    | 18    |
| MCM973D3.....  | 17    |
| MCM973D5.....  | 17    |

## MCO

|                 |    |
|-----------------|----|
| MCO201-2-5..... | 19 |
| MCO201-3-0..... | 19 |
| MCO305.....     | 20 |
| MCO770-1-5..... | 19 |
| MCO770-2-0..... | 19 |
| MCO770-2-5..... | 19 |
| MCO770-3-0..... | 19 |
| MCO770-3-5..... | 19 |
| MCO770-4-0..... | 19 |
| MCO771-2-5..... | 19 |
| MCO771-3-0..... | 19 |
| MCO771-4-0..... | 19 |
| MCO803.....     | 20 |
| MCO979.....     | 22 |
| MCO988.....     | 22 |

## MCR

|            |   |
|------------|---|
| MCR2A..... | 3 |
| MCR3A..... | 3 |

## P

|              |    |
|--------------|----|
| P761-7.....  | 50 |
| P761-8.....  | 50 |
| P762-2.....  | 51 |
| P762-4.....  | 51 |
| P763.....    | 52 |
| P764-1.....  | 51 |
| P764-3.....  | 51 |
| P764-5.....  | 51 |
| P764-7.....  | 51 |
| P764-11..... | 51 |
| P765-14..... | 52 |
| P765-16..... | 52 |
| P769.....    | 52 |
| P770.....    | 53 |
| P771.....    | 53 |

## S

|              |    |
|--------------|----|
| SM63(D)..... | 67 |
| SM64(D)..... | 67 |
| SM65(D)..... | 67 |

<sup>(D)</sup> Distributed by Integra MicroFrance SAS (France) - For corresponding manufacturer part numbers, see table in annex







**Table for corresponding  
manufacturer part numbers  
For information purposes only**

| Integra MicroFrance® product reference | Corresponding manufacturer part number |
|----------------------------------------|----------------------------------------|
| 1003IG                                 | M3IG                                   |
| 1003LI                                 | M3LI                                   |
| 1004IG                                 | M4IG                                   |
| 1005I                                  | M5I                                    |
| 7798T                                  | 101                                    |
| 7798T1                                 | 101 B                                  |
| 7872-24                                | 1224                                   |
| 7876-100                               | 310                                    |
| 7876-120                               | 312                                    |
| 7876-160                               | 316                                    |
| 7876-60                                | 306                                    |
| 7876-80                                | 308                                    |
| 7960T40                                | 731901                                 |
| MCL170 <sup>(1)</sup>                  | CO101                                  |
| MCL171 <sup>(1)</sup>                  | CO201                                  |
| MCL173 <sup>(1)</sup>                  | CO100                                  |
| MCL174 <sup>(1)</sup>                  | CO200                                  |
| MCL197 <sup>(1)</sup>                  | CO103                                  |
| MCL38N                                 | EC-BLU50U090                           |
| SM63                                   | SM63                                   |
| SM64                                   | SM64                                   |
| SM65                                   | SM65                                   |

<sup>(1)</sup> Not a medical device





Integra MicroFrance SAS (France)  
Le Pavillon ■ 03160 Saint Aubin le Monial ■ France



Swann Morton Ltd  
Penn Works ■ Owlerton Green ■ Sheffield  
South Yorkshire S6 2BJ ■ Royaume Uni



Emergo Europe  
Prinsessegracht 20  
2514 AP The Hague ■ The Netherlands

CE Class I products only :



Agencinox  
Zone Industrielle  
89110 Aillant sur Tholon ■ France



Directa AB  
Box 723  
19427 Upplands Vasby ■ Suède



Sunoptic Technologies  
6018 Bowdendale Avenue  
Jacksonville 32216 ■ Florida USA



AJW Technology Consulting GmbH  
Breite Straße 3  
40213 Düsseldorf ■ Germany

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Consult product labels and inserts for any indications, contraindications, hazards, warnings, precautions, and instructions for use.
- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.

#### Additional information for EMEA Customers only :

Products mentioned in this document are CE class I, Is, Ila or I Ib devices. Contact Integra should you need any additional information on devices classification.  
All the medical devices mentioned on this document are CE marked in accordance with the applicable European laws, unless specifically identified as "NOT CE MARKED".  
Products marked with the symbol<sup>(1)</sup> are distributed by Integra MicroFrance SAS (France)

Date of last revision : 03/2023

## For more information, please contact customer service :

### Europe, Middle-East, Africa

International: +33 (0) 437 47 59 50 ■ +33 (0) 437 47 59 25 (Fax) ■ [csemea@integralife.com](mailto:csemea@integralife.com)  
France: +33 (0) 437 47 59 10 ■ +33 (0) 437 47 59 29 (Fax) ■ [custservfrance@integralife.com](mailto:custservfrance@integralife.com)  
United Kingdom: +44 (0) 1264 312 725 ■ +44 (0) 1264 312 821 (Fax) ■ [custsvcs.uk@integralife.com](mailto:custsvcs.uk@integralife.com)  
Belgium & Luxembourg: +32 (0)2 257 4130 ■ +32 (0)2 253 2466 (Fax) ■ [custsvcbenelux@integralife.com](mailto:custsvcbenelux@integralife.com)  
Netherlands: +31(0) 852 083 167 ■ +31(0) 207 093 627 (Fax) ■ [custsvcnetherlands@integralife.com](mailto:custsvcnetherlands@integralife.com)  
Italy: +39 (0)2 577 8921 ■ +39 (0)2 575 113 71 (Fax) ■ [custsvcity@integralife.com](mailto:custsvcity@integralife.com)  
[integralife.eu/contact](http://integralife.eu/contact)

Integra and the Integra logo are registered trademarks of Integra LifeSciences Corporation or its subsidiaries in the United States and/or other countries.

MicroFrance is a trademark of Integra LifeSciences Corporation or its subsidiaries.

All other trademarks and product names are the property of their respective owners.

© 2023 Integra LifeSciences Corporation. All rights reserved.

0335699-5-EN

DOCUMENT FOR USE IN EUROPE, MIDDLE-EAST AND AFRICA ONLY.